


Measuring chart

mediven® for men

Standard size seamless

l=Lengths e=Circum-ferences



Lengths

| | I | II | III | IV | V | VI | VII |
|-----|-------|-------|-------|-------|-------|-------|-------|
| eD | 27-33 | 29-36 | 32-39 | 34-42 | 36-45 | 38-48 | 40-50 |
| eC | 28-34 | 30-37 | 33-40 | 35-43 | 37-46 | 39-49 | 41-51 |
| eB1 | 23-27 | 24-29 | 26-32 | 29-35 | 31-37 | 33-39 | 35-41 |
| eB | 18-20 | 20-22 | 22-24 | 24-26 | 26-28 | 28-30 | 30-32 |
| eY | 26-31 | 28-33 | 29-35 | 31-37 | 32-38 | 33-40 | 34-42 |
| eA | 17-22 | 19-24 | 21-26 | 23-29 | 25-32 | 27-34 | 28-35 |

Lengths

| | |
|-----------|----------|
| lD normal | 39-44 cm |
| lD petite | 34-38 cm |